**MANUFACTURER:**  
**ORTHOPEDIATRICS**  
2850 Frontier Drive  
Warsaw, IN 46582  
USA  
Ph 574-268-6379

## INSTRUMENT CARE, CLEANING & STERILIZATION INSTRUCTIONS

CE 0086

CI-0001 May 2013 Rev K, ECR 2509

### ENGLISH

### INSTRUMENT CARE, CLEANING & STERILIZATION INSTRUCTIONS

**Purpose/Scope:** This guide specifies the information for the processing of the medical devices claimed to be re-sterilizable and medical devices intended to be sterilized by the processor. This detailed information is provided so that the medical device can be processed safely and will continue to meet its performance specifications. Ref. ISO 17664:2004, ANSI/AAMI ST79:2006

**Device(s):** all reusable surgical instruments comprising fixed assemblies (no moving parts) and simple hinged assemblies and all single-use implantable devices supplied non-sterile, but intended to be used in a sterile state supplied by OrthoPediatrics Corp.

<b>WARNINGS</b>	<ol style="list-style-type: none"><li>Automated cleaning may not be effective. A thorough, manual cleaning process is recommended.</li><li>Where applicable, disassemble instruments prior to cleaning. Refer to Product IFU for instruments that require disassembly for cleaning.</li><li>Long narrow cannulations and blind holes require particular attention during cleaning.</li><li>Cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and must not be used. Enzymatic and cleaning agents with neutral pH are recommended.</li><li>Failure to properly clean the device may lead to inadequate sterilization.</li></ol>
<b>Limitations on Reprocessing</b>	Repeated processing, according to the instructions below, has minimal effect on OP reusable manual instruments or single use implants. End of life for instruments is normally determined by wear and damage due to use.

<b>INSTRUCTIONS</b>	
<b>Point of Use</b>	Remove excess body fluids and tissue with a disposable, non-shedding wipe and cover with a damp cloth. Body fluids and tissue should not be allowed to dry on instruments prior to cleaning.
<b>Containment/Transportation</b>	<ol style="list-style-type: none"><li>Universal precautions for handling contaminated/biohazardous materials should be observed.</li><li>Instruments should be cleaned within 30 minutes of use to minimize the potential for drying prior to cleaning.</li></ol>
<b>Preparation of Cleaning Agents</b>	Prepare neutral pH enzyme and cleaning agents at the use-dilution and temperature recommended by the manufacturer.
<b>Cleaning: Manual</b>	<ol style="list-style-type: none"><li>Use the neutral pH enzyme soaking solution that has been prepared.</li><li>Disassemble instruments, if applicable.</li></ol> <p>Completely submerge the instrument in enzyme solution and allow it to soak for a minimum of 10 minutes. Follow the enzymatic or cleaning agent manufacturer's instructions for the water quality and exposure time. Use a soft-bristled brush to gently clean the device (particular attention shall be given to crevices, lumens, mated surfaces and other hard-to-clean areas) until all visible soil has been removed. Lumens should be cleaned with a long, narrow, soft-bristled brush. Use a large syringe or pulsating water jet to thoroughly flush all channels and cavities with cleaning solution to remove debris. Repeat until all areas are clean.</p> <p><b>Note: The enzyme solution should be changed before it becomes grossly contaminated (bloody and/or turbid).</b></p> <ol style="list-style-type: none"><li>Remove the device from the enzyme solution and rinse in purified water (from one or any combination of the following processes: ultra-filter, RO, DI and/or distilled) for a minimum of 3 minutes. Thoroughly flush lumens, holes and other difficult to reach areas.</li><li>Visually inspect the device. Repeat the manual cleaning and rinsing procedure steps until no visible soil remains on the device.</li><li>Prepare the neutral pH cleaning (detergent) solution and place in a sonication unit.</li><li>Completely submerge device in cleaning solution and sonicate for 10 minutes, preferably at 45-50 KHz.</li><li>Rinse instrument in purified water (from one or any combination of the following processes: ultra-filter, RO, DI and/or distilled) thoroughly for at least 3 minutes or until there is no sign of blood or soil in the rinse stream.</li><li>Repeat step 5 with freshly prepared cleaning solution.</li><li>Repeat step 6 for thorough rinsing to remove any cleaning solution residues.</li><li>Dry the instrument with a clean, disposable, absorbent, non-shedding wipe or clean compressed air.</li></ol>
<b>Cleaning: Automated</b>	Automated washer/disinfecter systems are not recommended as the sole cleaning method for complex surgical instruments. These instruments should be cleaned following the manual cleaning procedure above. An automated system may be used as a follow-up method but is not required.
<b>Disinfection</b>	Disinfection is only acceptable as an adjunct to full sterilization for reusable surgical instruments. See sterilization section below.
<b>Drying</b>	Instruments must be thoroughly dried and all residual moisture must be removed before they are stored. Use a soft, absorbent towel/ cloth to dry external surfaces or clean compressed air.
<b>Inspection and Testing</b>	<ol style="list-style-type: none"><li>Carefully inspect each device to ensure that all visible blood and soil has been removed.</li><li>Visually inspect for damage and/or wear.</li><li>Check the action of moving parts (such as hinges and box-locks) to ensure smooth operation throughout the intended range of motion.</li><li>Check instruments with long slender features (particularly rotating instruments) for distortion.</li><li>Where instruments form part of a larger assembly, check that the devices assemble readily with mating components.</li></ol> <p><b>Note: Any instrument with heavy scratches, flaws, corrosion, cracked seals, discoloration, or which has been damaged or does not function properly has reached its end of life, and should be discarded following hospital protocol.</b></p>
<b>Maintenance</b>	Lubricate hinges, threads and other moving parts with a commercial water-based surgical grade instrument lubricant (such as instrument milk) to reduce friction and wear. Do not use mineral oil, petroleum, or silicone-based products. Discard blunt or damaged instruments.
<b>Packaging</b>	<ol style="list-style-type: none"><li>Singly – a standard polyethylene/Tyvek (or equivalent FDA cleared) sterilization pouch of the appropriate size may be used for single instruments. Ensure that the pack is large enough to contain the instrument without stressing the seals or tearing the packaging.</li><li>In Sets – sets of instruments may be loaded into dedicated instrument trays or general purpose sterilization trays for sterilization. If applicable, use standard medical grade steam sterilization FDA cleared wrap following the AAMI double wrap method (ANSI/AAMI ST79-2006).</li></ol>
<b>Sterilization</b>	<p><b>Warning: Failure to properly clean devices may lead to inadequate sterilization.</b></p> <p><b>OrthoPediatrics implants and instruments are recommended to be sterilized by steam autoclaving procedures regularly used in the hospital for wrapped instruments (based on ANSI/AAMI ST79: 2006) in accordance with the validated parameters. Sterilization cycle NOT validated for prion inactivation.</b></p> <p>Steam sterilize using validated method: Prevacuum cycle for 4 minutes at a temperature of 132°C Minimum drying time for Prevacuum cycle: 30minutes. Drying times will vary according to load size and should be increased for larger loads. The user should apply verifiable methods (e.g., visual inspection) to confirm adequate drying.</p> <p>The autoclave must be properly installed, maintained, and calibrated. Only legally marketed, FDA cleared sterilization wrap/pouches should be used by the end-user for packaging terminally sterilized devices. The autoclave manufacturer's operating instructions and recommended guidelines for maximum sterilization load should be followed.</p>
<b>Storage</b>	Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes.
<b>Additional Information</b>	1. When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer's maximum load is not exceeded. 2. Sterile instrument packages should be examined closely prior to opening to ensure that there has been no loss of package integrity.
<b>Customer Service Information</b>	OrthoPediatrics, Corp. 2850 Frontier Drive Warsaw, IN 46582 Phone: 574-268-6379 Fax: (574) 268-6302 Email: <a href="mailto:MarCom@OrthoPediatrics.com">MarCom@OrthoPediatrics.com</a>

It is the responsibility of the processor to ensure that reprocessing as actually performed using the appropriate equipment and materials, and personnel in the reprocessing facility have been adequately trained in order to achieve the desired result. This normally requires validation and routine monitoring of the process. Any deviation by the processor from these instructions should be properly evaluated for effectiveness to avoid potential adverse consequences. The recommended sterilization parameters are only valid with sterilization equipment that has been properly maintained.

This document is uncontrolled after distribution and it is the responsibility of the reader to confirm or obtain the most current revision of this document by contacting OrthoPediatrics Customer Service

OrthoPediatrics and the Pedi logo are trademarks, registered in the United States, of OrthoPediatrics Corp.

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OrthoPediatrics et le logo Pedi sont des marques de commerce déposées aux États-Unis d'OrthoPediatrics Corp.
Logo Pedi dan OrthoPediatrics adalah tanda dagangan OrthoPediatrics Corp., didaftarkan di Amerika Syarikat.
OrthoPediatrics en het Pedi-logo zijn in de Verenigde Staten gedeponeerde handelsmerken van OrthoPediatrics Corp.